



# INSTRUCTIONS FOR USE

## Oosafe® Hand Disinfectant

EN-BIO-OODH-01.05.2024

# DISINFECTANTS



**Oosafe® Hand Disinfectant**  
 REF OODH-0500

**ME and HSSA tested for each LOT**

Hand disinfectant solution, non-sensitizing, non-irritating and suitable for frequent use.  
 Contains alcohol, odorless and does not release volatile organic compounds.  
 Bactericidal and Virucidal on Vaccinia virus (all enveloped viruses from Annex A, EN14476:2013+A1:2015 including HBV, HCV, HDV, HIV, Coronavirus, Influenza, Rabies, Measles virus).  
 Fungicidal and sporicidal on Bacillus spores.  
 Follow laboratory's hygienic procedure for washing hands before using hand disinfectant. Apply minimum 2 pumps (3 ml) of disinfectant and rub hands until they are completely dry.  
 Suitable for aquatic life with long-lasting effects. Avoid release to the environment. Dispose of contents in accordance with local regulations. Keep out of the reach of children.  
 Contains: Quaternary ammonium compounds; Benzyl-alkylmethyl chloride 2.6 mg/g.

LOT 81010 ID 87303(A)

**500 ml**

SparMED ApS  
 Ryttersmarken 2,  
 4291 SØ,  
 Denmark  
 www.sparmed.dk

5 712502 400046

**Oosafe® Hand Disinfectant**  
 REF OODH-5000

**5 L**

2024-05 LOT 81000 ID 803(5A)

SparMED ApS  
 Ryttersmarken 2,  
 4291 SØ,  
 Denmark  
 www.sparmed.dk

5 712502 400046



EN-BIO-ODDH-01.05.2024

## Attention!

- Check the expiration date on the label;
- Storage temperature is between 15 – 30°C;
- Avoid direct sunlight;
- Visually check for bottles to be undamaged and the sealing intact;
- Harmful to aquatic life with long-lasting effects;
- Avoid release to the environment;
- Keep out of the reach of children;
- Oosafe® disinfectant bottles are made of high-density polyethylene (HDPE) and are recyclable. For the sake of the environment, it is recommended to dispose bottles according to the local regulations.



- Dispose of contents in accordance with local regulations.
- Cardboard boxes shall be disposed according to the local regulations for paper recycling;
- Plastic bottles and cardboard boxes are flammable. Keep away from fire.




## First-aid measures!

- Inhalation: Supply fresh air; consult doctor in case of complaints.
- Skin contact: The product is non-sensitizing and non-irritating, suitable for frequent use.
- Eye contact: Rinse opened eyes for 10 minutes under running water. Force to open eyelids. Consult a doctor if irritation persists.
- Ingestion: Rinse mouth and drink plenty of water. Do not cause vomiting. If symptoms persist, consult doctor.

## Additional information:

For more detailed information about the products, safety refer to disinfectants' safety data sheet.

**DISIN  
FECT  
ANTS**



### Intended use

- Hand disinfection.

### Product description

Oosafe® Hand Disinfectant is a colourless, odourless, ready to use water-based hand disinfectant. The disinfectant is alcohol-free, non-fragranced, does not release volatile organic compounds (VOCs) and therefore is a favourable disinfectant for use by professionals in in-vitro fertilization (IVF) laboratories.

### Application

Follow laboratory's hygienic procedure for washing hands before using hand disinfectant.

Use a clean knife or scissor to open the seal and remove it completely.


For effective hand disinfection, apply minimum 3 ml of disinfectant and rub hands until they are completely dry.

### Features

- Clear, odourless;
- Ready-to-use;
- Alcohol-free;
- Non-VOC releasing;
- Mouse embryo assay (MEA) and Human Sperm Survival Assay (HSSA) tested for each LOT;
- Bactericidal, Yeastocidal and Virucidal on Vaccinia virus.

### Durability

After opening Oosafe® Disinfectant bottles or containers, the mixture maintains stability for 3 years under proper storage condition.



## Microbiological efficacy

Efficacy table for Hand Disinfectant		
Activity Spectrum	Method*	Contact Time
<b>Bactericidal</b> <i>Pseudomonas aeruginosa,</i> <i>Staphylococcus aureus, Enterococcus hirae,</i> <i>Escherichia coli</i>	EN 1040 EN 13727+A2 EN 14561 EN 1500	30 sec
<b>Yeasticidal</b> <i>Candida albicans</i>	EN 1275 EN 13624 EN 14562	30 sec
<b>Virucidal</b> <i>Vaccinia virus**</i>	EN 14476+A1	30 sec

\*Efficacy is tested for clean conditions.

\*\*Effective against all enveloped viruses from Annex A, EN14476:2013+A1:2015 including HBV, HCV, HDV, HIV, Coronavirus, Influenza, Rabies, Measles virus.









## Packaging Type

Packaging Type	Product reference
1. 50 ml, spray bottle - sample	OODH-0050
2. 500 ml, bottle with pump: individually packed in the carton box	OODH-0500
3. 5 L, refill container, individually packed in the cardboard protection	OODH-5000

Symbols on the product's primary packing label

**DISIN  
FECT  
ANTS**



Symbol	Description of the symbol
1. 	Consult instructions for use or consult electronic instructions for use
2. 	Manufacturer
3. 	Use by date
4. 	Temperature limit
5. 	Batch code
6. 	Keep away from sunlight
7. 	Identification number
8. 	Catalogue number (Product reference)

Manufacturer:

**sparomed**

**SparMED ApS**

Ryttermarken 2, DK-3520 Farum, Denmark

Phone: +45 3940 2503, E-mail:

[info@sparmed.dk](mailto:info@sparmed.dk) Website: [www.sparmed.dk](http://www.sparmed.dk)

