

Oosafe® AIR CLEANER KN



CE Manufactured
in accordance with
ISO 9001:2015

Oosafe® AIR CLEANER KN

Oosafe® Air Cleaner KN is advanced equipment, designed for IVF laboratory use. It purifies air against volatile organic compounds (VOCs), chemically active compounds (CACs), microorganisms and particles.



- **Better air quality:** Oosafe® Air Cleaner KN purifies the air in laboratories.
- **More activated carbon, more VOC removal capacity:** Canister Carbon Filter Set consists of 5 canisters made of pure activated carbon.
- **45° air blowing:** Faster and more efficient air purification with 45° air blowing system.
- **Low noise level:** Higher user comfort is achieved with very low noise level (lower than 50dB (A) at full speed and 42dB (A) at half speed).
- **Known worldwide:** Oosafe® Air Cleaner KN has been used for many years in well-known IVF clinics all over the world as the main air purification device.



Oosafe® Air Cleaner KN

Product code:
OOACEQ-03

Oosafe® AIR CLEANER KN

OOSAFE® AIR CLEANER KN

- **Adjustable and Accurate Air Speed Level:**

Air speed can be adjusted from the LCD Screen.

- **Alarm functions:**

Timer alerts when filters should be replaced. No replacement of filters before it is needed.

- **Easy to operate:**

Main service can be handled on the LCD screen through service menu. Replacement of the filters can be done by the users.

- **Large HEPA filter:**

Does not need to be replaced as often. Better air quality is achieved by purifying airborne particles and microorganisms.

TECHNICAL DATA

Dimensions	660x533x1640 (w x d x h, in mm)
Air exchanges per hour (AEH)	9-28 times (varies on size of the room)
Weight (including filters)	77 kg
HEPA Filter efficiency (MPPS)	≥ 99.995%
Total amount of activated carbon	9.4 kg

Oosafe® Air Cleaner KN

Product code:
OOACEQ-03



ORDER INFORMATION


Product code	Product name
OOACEQ-03	Oosafe® Air Cleaner KN without initial filter set
OOACHF-02	HEPA Filter for Oosafe® Air Cleaner KN
OOACCF-01	Canister Carbon Filter Set for Oosafe® Air Cleaner KN
OOACPF-01	Main Prefilter for Oosafe® Air Cleaner KN

Oosafe® AIR CLEANER KN



HEPA Filter
for Oosafe® Air Cleaner KN

Product code:
OOACHF-02



Canister Carbon Filter Set
for Oosafe® Air Cleaner KN

Product code:
OOACCF-01



Main Prefilter
for Oosafe® Air Cleaner KN

Product code:
OOACPF-01

sparomed



SparMED ApS
Ryttermarken 2, 3520 Farum, DENMARK
Phone number: +45 39-40-25-03



MK_4.3.B8 rev. 20.06.2023

